

AED Inspection Log

School/Campus: _____

Month/Year: _____

AED unit/ID: _____

Check each AED regularly (many schools inspect monthly). Keep this log on file as part of your CERP documentation.

Date	AED location	Ready light OK	Pads exp.	Battery exp.	Damage? (Y/N)	Inspector

NOT SURE WHERE YOUR CAMPUS STANDS?

Book a free School Rescue Ready™ Review. We walk your campus, time your AED coverage against the three-minute standard, and show you exactly what is missing. No cost, no obligation, and you keep what we find either way.

Call or text (919) 372-9657 · info@rescuebeatspro.com